



mustangTM

The Mustang is an anterior/posterior gait trainer that provides proper positioning and support for children while they learn the skills of stepping/walking



The Mustang provides excellent upper body support and adjustable angle positioning for forward tilt in the anterior configuration and full upright support in the posterior configuration. A variety of positioning accessories are available for children who need extra support, and are removable as their walking skills progress. The center spar is dynamic and helps mimic the natural up and down movement in walking.

The Mustang is an excellent choice for children with Cerebral Palsy and other developmental delays. It offers flexibility, growth and ease of use while children are learning to walk and adapting to their changing needs.

The Mustang comes in four sizes suitable for children from 1 to 18 years of age.



Motivating the child to walk

The user is placed at the center of the Mustang which provides stability of the upper body and enables the user to focus energy on walking



Develop a more natural walking pattern

The Mustang has an easy and steady angle adjustable center bar. A slight forward leaning position, combined with built-in suspension, stimulates the natural walking pattern



Easy to access

The Mustang has been constructed with accessibility in mind. The U-shaped frame and wrap around support, allows for easy access and transfer



Easy to store and carry

The foldable frame makes storage and transport of the Mustang easy



Standard

As standard the Mustang is delivered with lockable wheels. The Mustang is available in 4 sizes



Size 1



Size 2, 3 and 4 with larger wheels

Frame color

Available in 4 colors



Size 1 = Yellow



Size 2 = Blue



Size 3 = Red



Size 4 = Blue

Wheels

The Mustang is delivered as standard with various adjustment possibilities on all four wheels



Directional stabilizer



Tension adjustable wheels



Wheel brake on all 4 wheels




Lock



Release

Wheel stops

The unique on/off anti-reverse wheel stops are very easy to lock and unlock



Measurements

	Size 1	Size 2
Width, total (G)	52 cm (20 ¹ / ₄ ")	63 cm (24 ¹ / ₂ ")
Width, inside frame (G1)	38 cm (14 ³ / ₄ ")	45 cm (17 ¹ / ₂ ")
Length, total (H)	75 cm (29 ¹ / ₄ ")	90 cm (35")
Height, total (I)	65 cm (25 ¹ / ₄ ")	80 cm (31 ¹ / ₄ ")
Seat height above floor (seat = accessory)	17-35 cm (6 ¹ / ₂ -13 ¹ / ₂ ")	33-50 cm (12 ³ / ₄ -19 ¹ / ₂ ")
Width between handles (handles = accessory)	14-35 cm (5 ¹ / ₂ -13 ¹ / ₂ ")	14-35 cm (5 ¹ / ₂ -13 ¹ / ₂ ")
Height from floor to handles (handles = accessory)	37-58 cm (14 ¹ / ₂ -22 ¹ / ₂ ")	51-71 cm (20-27 ³ / ₄ ")
Angle range of center bar (a1)	30°	30°
Turning radius	0°	0°
User height, floor to axil	38-65 cm (14 ³ / ₄ -25 ¹ / ₄ ")	60-80 cm (23 ¹ / ₂ -31 ¹ / ₄ ")
Weight, standard without accessories	10 kg (22 lb)	14 kg (31 lb)
Max. load / user weight	30 kg (66 lb)	40 kg (88 lb)
	Size 3	Size 4
Width, total (G)	65 cm (25 ¹ / ₂ ")	70 cm (27 ¹ / ₂ ")
Width, inside frame (G1)	50 cm (19 ¹ / ₂ ")	56 cm (22")
Length, total (H)	100 cm (39")	110 cm (43")
Height, total (I)	105 cm (41")	130 cm (50 ³ / ₄ ")
Seat height above floor (seat = accessory)	45-68 cm (17 ¹ / ₂ -26 ¹ / ₂ ")	62-84 cm (24-32 ³ / ₄ ")
Width between handles (handles = accessory)	14-60 cm (5 ¹ / ₂ -23 ¹ / ₂ ")	14-60 cm (5 ¹ / ₂ -23 ¹ / ₂ ")
Height from floor to handles (handles = accessory)	65-90 cm (25 ¹ / ₂ -35")	81-125 cm (31 ¹ / ₂ -48 ³ / ₄ ")
Angle range of center bar (a1)	0-25°	0-25°
Turning radius	55 cm (21 ¹ / ₂ ")	62 cm (24 ³ / ₈ ")
User height, floor to axil	70-95 cm (27 ¹ / ₂ -37")	87-130 cm (34-50 ³ / ₄ ")
Weight, standard without accessories	18 kg (40 lb)	19 kg (42 lb)
Max. load / user weight	60 kg (132 lb)	80 kg (176 lb)

The Mustang is measured as standard delivered and with the center bar in vertical position. Different settings may result in another set of measurements. To accomplish some of the measurements above, a slight adjustment might be required. We supply various types of accessories, which is why the weight is listed without.



ACCESSORIES



Chest support

8698x0

For support of the upper body

size	measurements
1	O: 15½"-25¼"
2	O: 26¼"-30¾"
3	O: 34¼"-44¾"

1

O: 15½"-25¼"

2

O: 26¼"-30¾"

3

O: 34¼"-44¾"

O: inner circumference



Cover

8698x0-x60

For added comfort

size	measurements
1	H:3¼" x W:9¾"
2	H:4¾" x W:14¾"
3	H:5¾" x W:19"

1

H:3¼" x W:9¾"

2

H:4¾" x W:14¾"

3

H:5¾" x W:19"



Extension piece

869225

To eliminate the gap when width adjusting the chest support

size	measurements
one	H:4¾" x W:2¼"

one

H:4¾" x W:2¼"



Hand bar

86936x

Provide hand control For anterior use

size	measurements
S	H:4¼" x W:8¼"
L	H:4¼" x W:5½"

S

H:4¼" x W:8¼"

L

H:4¼" x W:5½"



Abdominal support

869260

Provides additional anterior pelvic control

size	measurements
one	H:5" x W:7½"

one

H:5" x W:7½"



Arm rests, anterior

86935x. Provide forearm support. Helps the child to stabilize the trunk and head

Two sizes



Sacrum support

86925x

Helps to stabilize the pelvis. Requires seat 8691X0

size	measurements
1	H:7" (from seat) D:5-8¼" (from seat)
2	H:8" (from seat) D:6½"-9½" (from seat)

1

H:7" (from seat)
D:5-8¼" (from seat)

2

H:8" (from seat)
D:6½"-9½" (from seat)

Arm prompts, anterior or posterior

8693xx

Helps the child to stabilize the trunk and head

Three sizes



Hip support, small

869220 or 869221

Helps to stabilize the pelvis. Requires seat 8691X0

size	measurements
one	W:7¼"-10"

one

W:7¼"-10"



Hip support, large

869230 or 869231

Helps to stabilize the pelvis. Requires seat 8691X0

size	measurements
one	W:10-12¾"

one

W:10-12¾"



Head support system

861424-3 and 861424-4

size	measurements
1	Rod = L:7¾"
2	Rod = L:12¾"

1

Rod = L:7¾"

2

Rod = L:12¾"



Seat

8691x0. Provides a resting position that a child can use when fatigued. Used for mounting sacrum or hip supports

size	measurements
1	H:6¾" x D:10"
2	H:7" x D:12¾"

1

H:6¾" x D:10"

2

H:7" x D:12¾"



Connecting rod

869226

Mounted between hip and head support

One size



Sling seat

86957x

Provides partial weight bearing support

size	measurements
1	W:7 x L:9½"
2	W:7 x L:12¾"



Leg pad

869470

Provides knee protection

One size



Guide bar

869510

Allows caregiver to direct and assist the child's movements

size	measurements
one	L: 27¼"



Storage pouch

869560

For transporting books etc.
Only for size 1+2

size	measurements
one	H:5¼ x W:5"



Storage bag

869565

For transporting books etc.
Only for size 3+4

size	measurements
one	H:12 x W: 15½"



Leg separator

8694x0

Protects children with the tendency of scissoring their legs

Four sizes



Ankle prompts

8694xx

Guidance for correct walking pattern

size	measurements
1	L: 20½"
2	L: 24 & 10½"
3	L: 27¾"
4	L: 31¾"